






































































	IUPAC Name	Synonyms		Linear formula	CAS №	Molecular weight	Assay %	Density g/cm <sup>3</sup>	Refractive index n <sub>20</sub> /D	Melting point °C	Boiling point °C
		Синонимы	Наименование	Структурная формула	CAS №	Молекулярный вес	Содержание основного компонента %	Плотность g/cm <sup>3</sup>	Коэффициент преломления n <sub>20</sub> /D	Температура плавления °C	Температура кипения °C
	1,1,2,3,4,4-Hexafluoro-1,3-butadiene	Hexafluoro-1,3- butadiene	1,1,2,3,4,4-Гексафтор-1,3-бутадиен	C <sub>4</sub> F <sub>6</sub>	685-63-2	162,033	99,995	1,553 (T= -20°C)		-132,1	6-7
	Trimethylsilane	3MS	Триметилсилан	C <sub>3</sub> H <sub>10</sub> Si	993-07-7	73,18996	99,999	0,634 (T= 6,7 °C)		-135,9	6,7
	Bromo(trimethyl)silan	Trimethylbromosilane Me <sub>3</sub> SiBr Trimethylsilyl bromide	Бром(триметил)силан	C <sub>3</sub> H <sub>9</sub> BrSi	2857-97-8	153,09		1,16 (T=25 °C)	1,424 (lit)	-43	79 (lit)
	Hexafluorobenzene	Perfluorobenzene	Гексафторбензол	C <sub>6</sub> F <sub>6</sub>	392-56-3	186,05	99,99	1,612 (T=25 °C)	1,377	3,7-4,1 (lit)	80-82 (lit)
	Bromoform	Tribromomethane	Бромформ	CHBr <sub>3</sub>	75-25-2	252,73094	98	2,89		8	150-151
	Tetrakis (dihydrocarbylamino) phosphonium bromide			C <sub>16</sub> H <sub>40</sub> BrN <sub>4</sub> P	81175-49-7						
	Aluminium tribromide	Aluminium bromide	Бромид алюминия	AlBr <sub>3</sub>	7727-15-3	266,69354	99	3,2		97	265
	Tetrachlorostannane	Tin tetrachloride anhydrous Tin (IV) chloride, anhydrous	Олово четыреххлористое безводное	SnCl <sub>4</sub>	7646-78-8	260,502	99	2,22		-33 (lit)	114 (lit)
	Tetrachlorostannane pentahydrate	Tin tetrachloride pentahydrate Tin (IV) chloride pentahydrate	Олово четыреххлористое пятиводное	SnCl <sub>4</sub> ·5H <sub>2</sub> O	10026-06-9	350,5784	99,9			56	
	1,2,3,4,5-Pentafluoro-6-(trifluoronethyl)benzene	Octafluorotoluene	Октафтортолуол	C <sub>6</sub> F <sub>5</sub> CF <sub>3</sub>	434-64-0	236,06	99	1,666 (T=25 °C)	1,368 (lit)		104 (lit)

	1,2,3,4-Tetrafluoro-5-(trifluoromethyl)benzene	2,3,4,5-Tetrafluorobenzotrifluoride	1,2,3,4-Тетрафтор-5-(трифторметил)бензол	C <sub>7</sub> HF <sub>7</sub>	654-53-5	218,073761	98				104
	1,2,4,5-Tetrafluoro-3-(trifluoromethyl)benzene	2,3,5,6-Tetrafluorobenzotrifluoride	1,2,4,5-Тетрафтор-3-(трифторметил)бензол	C <sub>7</sub> HF <sub>7</sub>	651-80-9	218,073761	99	1,601 (T=20°C)	1,379		111-112
	1,2,3,5-Tetrafluoro-4-(trifluoromethyl)benzene	2,3,4,6-Tetrafluorobenzotrifluoride	1,2,3,5-Тетрафтор-4-(трифторметил)бензол	C <sub>7</sub> HF <sub>7</sub>	5360-82-7	384,53652	99				
	1,2,3-Trifluoro-4-(trifluoromethyl)benzene	2,3,4-Trifluorobenzotrifluoride	1,2,3-Трифтор-4-(трифторметил)бензол	C <sub>7</sub> H <sub>2</sub> F <sub>6</sub>	393-01-1	200,083298	98			1,383-1,385	
	1,3,5-Trifluoro-2-(trifluoromethyl)benzene	2,4,6-Trifluorobenzotrifluoride	1,3,5-Трифтор-2-(трифторметил)бензол	C <sub>7</sub> H <sub>2</sub> F <sub>6</sub>	122030-04-0	200,083298	99	1,475 (T=20°C)	1,384-1,386		105-107
	1,2,3-Trifluoro-5-(trifluoromethyl)benzene	3,4,5-Trifluorobenzotrifluoride	1,2,3-Трифтор-5-(трифторметил)бензол	C <sub>7</sub> H <sub>2</sub> F <sub>6</sub>	80172-04-9	200,083298	98				
	Pentafluorobenzoic acid	Perfluorobenzoic acid	Пентафторбензойная кислота	C <sub>7</sub> HF <sub>5</sub> O <sub>2</sub>	602-94-8	212,07	99	1,944 (T=25°C)		100-102 (lit)	220 (lit)
	2,3,4,5-Tetrafluorobenzoic acid		2,3,4,5-Тетрафторбензойная кислота	C <sub>7</sub> H <sub>2</sub> F <sub>4</sub> O <sub>2</sub>	1201-31-6	194,085	98			85-87	
	2,3,5,6-Tetrafluorobenzoic acid		2,3,5,6-Тетрафторбензойная кислота	C <sub>7</sub> H <sub>2</sub> F <sub>4</sub> O <sub>2</sub>	652-18-6	194,085	99			148-152	
	2,3,4,6-Tetrafluorobenzoic acid		2,3,4,6-Тетрафторбензойная кислота	C <sub>7</sub> H <sub>2</sub> F <sub>4</sub> O <sub>2</sub>		194,085	98				
	2,3,4-Trifluorobenzoic acid		2,3,4-Трифторбензойная кислота	C <sub>7</sub> H <sub>3</sub> F <sub>3</sub> O <sub>2</sub>	61079-72-9	176,09	98			140-142	
	2,4,6-Trifluorobenzoic acid		2,4,6-Трифторбензойная кислота	C <sub>7</sub> H <sub>3</sub> F <sub>3</sub> O <sub>2</sub>	28314-80-9	176,09	99			138-144	
	3,4,5-Trifluorobenzoic acid		3,4,5-Трифторбензойная кислота	C <sub>7</sub> H <sub>3</sub> F <sub>3</sub> O <sub>2</sub>	121602-93-5	176,09	98			95-99	
	4-Chloro-2,3,5,6-tetrafluorobenzoic acid		4-Хлор-2,3,5,6-тетрафторбензойная кислота	C <sub>7</sub> HCIF <sub>4</sub> O <sub>2</sub>		228,528	98	1,7	1,494		257,4±35
	4-Bromo-2,3,5,6-tetrafluorobenzoic acid	2,3,5,6-tetrafluoro-4-bromobenzoic acid	4-Бром-2,3,5,6-тетрафторбензойная кислота	C <sub>7</sub> HBrF <sub>4</sub> O <sub>2</sub>	4707-24-8	271,981352	98			144-145	
	2-Bromo-3,4,5,6-tetrafluorobenzoic acid		2-Бром-3,4,5,6-тетрафторбензойная кислота	C <sub>7</sub> HBrF <sub>4</sub> O <sub>2</sub>	16583-04-3	271,981352	98	2,0 ± 0,1	1,517		266,7±40
	3-Bromo-2,4,5,6-tetrafluorobenzoic acid		3-Бром-2,4,5,6-тетрафторбензойная кислота	C <sub>7</sub> HBrF <sub>4</sub> O <sub>2</sub>	33564-64-6	271,981352	98	2,0 ± 0,1	1,517		262,9±40

	2,3,5,6-Tetrafluoro-4-methylbenzoic acid	4-Methyl-2,3,5,6-tetrafluorobenzoic acid 2,3,5,6-Tetrafluoro-p-toluic acid	2,3,5,6-тетрафтор- 4-метилбензойная кислота	C <sub>8</sub> H <sub>4</sub> F <sub>4</sub> O <sub>2</sub>	652-32-4	208,112172				168-172	
	(Pentafluorophenyl)methanol	Pentafluorobenzyl alcohol	Пентафторфенил метиловый спирт	C <sub>7</sub> H <sub>3</sub> F <sub>5</sub> O	440-60-8	198,092235	98			37-38	114-115 (p=60 torr)
	(2,3,4,5-Tetrafluorophenyl)methanol	2,3,4,5-Tetrafluorobenzyl alcohol	(2,3,4,5-Тetraфторфенил) метиловый спирт	C <sub>7</sub> H <sub>4</sub> F <sub>4</sub> O	53072-18-7	180,1011772	98		1,454-1,456		140 (p=90 torr)
	(2,3,5,6-Tetrafluorophenyl)methanol	2,3,5,6-Tetrafluorobenzyl alcohol	(2,3,5,6-Тetraфторфенил) метиловый спирт	C <sub>7</sub> H <sub>4</sub> F <sub>4</sub> O	4084-38-2	180,1011772	98			32-38	220-225
	(2,3,5,6-Tetrafluoro-4-methylphenyl)methanol	4-Methyl-2,3,5,6-tetrafluorobenzyl alcohol	(2,3,5,6-Тetraфтор-4-метилфенил) метиловый спирт	C <sub>8</sub> H <sub>6</sub> F <sub>4</sub> O	79538-03-7	194,128652				48	
	(2,3,4,6-Tetrafluorophenyl)methanol	2,3,4,6-Tetrafluorobenzyl alcohol	(2,3,4,6-Тetraфторфенил) метиловый спирт	C <sub>7</sub> H <sub>4</sub> F <sub>4</sub> O	53001-70-0	180,100	98	1,5 ±0,1	1,457		180,2±35
	(2,3,4-Trifluorophenyl)methanol	2,3,4-Trifluorobenzyl alcohol	(2,3,4-Трифторфенил) метиловый спирт	C <sub>7</sub> H <sub>5</sub> F <sub>3</sub> O	144284-24-2	162,111309			1,47	29,6	
	(2,4,6-Trifluorophenyl)methanol	2,4,6-Trifluorobenzyl alcohol	(2,4,6-Трифторфенил) метиловый спирт	C <sub>7</sub> H <sub>5</sub> F <sub>3</sub> O	118289-07-9	162,111309				91-94	230
	(3,4,5-Trifluorophenyl)methanol	3,4,5-Trifluorobenzyl alcohol	(3,4,5-Трифторфенил) метиловый спирт	C <sub>7</sub> H <sub>5</sub> F <sub>3</sub> O	220227-37-2	162,111309			1,4725		203-204
	(3,4-Difluorophenyl)methanol	3,4-Difluorobenzyl alcohol	(3,4-Дифторфенил) метиловый спирт	C <sub>7</sub> H <sub>6</sub> F <sub>2</sub> O	85118-05-4	144,120846		1,26	1,4885-1,4905		200
	1-(2,3,4,5-Tetrafluorophenyl)methamine	2,3,4,5-Tetrafluorobenzyl amine	1-(2,3,4,5-Тetraфторфенил) метиламин	C <sub>7</sub> H <sub>5</sub> F <sub>4</sub> N	235088-70-7	179,115					
	1-(2,3,5,6-Tetrafluorophenyl)methanamine	2,3,5,6-Tetrafluorobenzyl amine	1-(2,3,5,6-Тetraфторфенил)метиламин	C <sub>7</sub> H <sub>5</sub> F <sub>4</sub> N	89992-52-9	179,115					
	1-(2,3,4,6-Tetrafluorophenyl)methamine	2,3,4,6-Tetrafluorobenzyl amine	1-(2,3,4,6-Тetraфторфенил) метиламин	C <sub>7</sub> H <sub>5</sub> F <sub>4</sub> N		179,115	98				
	1-(2,4,6-Trifluorophenyl)methanamine	2,4,6-Trifluorobenzyl amine	1-(2,4,6-Трифторфенил)метиламин	C <sub>7</sub> H <sub>6</sub> F <sub>3</sub> N	214759-21-4	161,126549					165
	1-(2,4,5-Trifluorophenyl)methanamine	3,4,5-Trifluorobenzyl amine	1-(2,4,5-Трифторфенил)метиламин	C <sub>7</sub> H <sub>6</sub> F <sub>3</sub> N	235088-69-4	161,126549			1,4736		

	1-(3,4-Difluorophenyl)methylamine	3,4-Difluorobenzyl amine	1-(3,4-Дифторфенил)метиламин	C <sub>7</sub> H <sub>7</sub> F <sub>2</sub> N	72235-53-1	143,136086		1,21	1,492-1,494		
	Pentafluorophenol		Пентафтофенол	C <sub>6</sub> HF <sub>5</sub> O	771-61-9	184,065355	99	1,757	1,427	34-36	143
	2,3,5,6-Tetrafluorophenol		2,3,5,6-Тетрафторфенол	C <sub>6</sub> H <sub>2</sub> F <sub>4</sub> O	769-39-1	166,074892	99			33-38	140
	Pentafluoroaniline	2,3,4,5,6-Pentafluorobenzeneamine	Пентафторанилин	C <sub>6</sub> H <sub>2</sub> F <sub>5</sub> N	771-60-8	183,080595	98	1,744 (T=20 °C)		33-38	153
	2,3,5,6-Tetrafluoroaniline	1-Amino-2,3,5,6-tetrafluorobenzene	2,3,5,6-Тетрафторанилин	C <sub>6</sub> H <sub>3</sub> F <sub>4</sub> N	700-17-4	165,090132		1,744 (T=20 °C)		29-33	158
	2,3,4,6-Tetrafluoroaniline	2,3,4,6-Tetrafluorobenzeneamine	2,3,4,6-Тетрафторанилин	C <sub>6</sub> H <sub>3</sub> F <sub>4</sub> N	363-73-5	165,090132		1,744 (T=20 °C)		23-26	65 (p=20 torr)
	2,3,4-Trifluoroaniline	2,3,4-Trifluorobenzeneamine	2,3,4-Трифторанилин	C <sub>6</sub> H <sub>4</sub> F <sub>3</sub> N	3862-73-5	147,099669		1,39	1,487		92 (p=48 torr)
	2,4,6-Trifluoroaniline	2,4,6-Trifluorobenzeneamine	2,4,6-Трифторанилин	C <sub>6</sub> H <sub>4</sub> F <sub>3</sub> N	363-81-5	147,099669		1,39		33-37	57 (p=22 torr)
	2,3,6-Trifluoroaniline	2,3,6-Trifluorobenzeneamine	2,3,6-Трифторанилин	C <sub>6</sub> H <sub>4</sub> F <sub>3</sub> N	67815-56-9	147,099669		1,39	1,486-1,488		150-151
	2,3,5-Trifluoro-6-nitroanilin	3,5,6-Trifluoro-2-nitroaniline	2,3,5-Трифтор-6-нитроанилин	C <sub>6</sub> H <sub>3</sub> F <sub>3</sub> N <sub>2</sub> O <sub>2</sub>	5415-62-3	192,097229					
	3,4,6-Trifluor-1,2-benzendiamine	1,2-Diamino-3,4,6-trifluorobenzene	3,4,6-Трифтор-1,2-бензолдиамин	C <sub>6</sub> H <sub>5</sub> F <sub>3</sub> N <sub>2</sub>		162,113					
	1,2,3,5-Tetrafluorobenzene	1,2,4,6-Tetrafluorobenzene	1,2,3,5-Тетрафторбензол	C <sub>6</sub> H <sub>2</sub> F <sub>4</sub>	2367-82-0	150,075492		1,31	1,4041	-48	83
	1,2,4,5-Tetrafluorobenzene	2,3,5,6-Tetrafluorobenzene	1,2,4,5-Тетрафторбензол	C <sub>6</sub> H <sub>2</sub> F <sub>4</sub>	327-54-8	150,075492		1,42	1,407	4	90
	1,2,3,4,5-Pentafluorobenzene	Pentafluorobenzene	1,2,3,4,5-Пентафторбензол	C <sub>6</sub> HF <sub>5</sub>	363-72-4	168,065955		1,51	1,39	-48	85
	1-Brom-2,3,4,5,6-pentafluorobenzene	Bromopentafluorobenzene	1-Бром-2,3,4,5,6-пентафторбензол	C <sub>6</sub> BrF <sub>5</sub>	344-04-7	246,962015		1,981	1,449	-31	137
	2,3,4,5,6-Pentafluoropyridine	Pentafluoropyridine	2,3,4,5,6-Пентафторпиридин	C <sub>5</sub> F <sub>5</sub> N	700-16-3	169,053715	99	1,54	1,386-1,387	-41,5	84

	2,3,5,6-Tetrafluoropyridine		2,3,5,6-Тетрафторпиридин	C <sub>5</sub> HF <sub>4</sub> N	2875-18-5	151,063252		1,52	1,4046	102	102
	2,3,6-Trifluoropyridine		2,3,6-Трифторпиридин	C <sub>5</sub> H <sub>2</sub> F <sub>3</sub> N	3512-18-3	133,072789		1,499	1,42		100-102
	2,3,4,6-Tetrafluoropyridine		2,3,4,6-Тетрафторпиридин	C <sub>5</sub> HF <sub>4</sub> N	3512-13-8	151,063252			1,399		88-90
	2,6-Difluoropyridine		2,6-Дифторпиридин	C <sub>5</sub> H <sub>3</sub> F <sub>2</sub> N	1513-65-1	115,082326		1,26	1,436-1,438		124,5 (p=743 torr)
	3-Chloro-2,4,5,6-tetrafluoropyridine	3-Chlorotetrafluoropyridine	3-Хлор-2,4,5,6-тетрафторпиридин	C <sub>5</sub> ClF <sub>4</sub> N	1735-84-8	185,508312		1,64	1,4345-1,4365		118-121
	3-Chloro-2,5,6-trifluoropyridine		3-Хлор-2,5,6-трифторпиридин	C <sub>5</sub> HClF <sub>3</sub> N	2879-42-7	167,517849			1,456		132-134
	2,4,6-Trifluoropyridine		2,4,6-Трифторпиридин	C <sub>5</sub> H <sub>2</sub> F <sub>3</sub>	3512-17-2	133,072789		1,499 (T=20 °C)	1,412		96,5-102
	2,3,5-Trifluoropyridine		2,3,5-Трифторпиридин	C <sub>5</sub> H <sub>2</sub> F <sub>3</sub>	76469-41-5	133,072789	98	1,499 (T=20 °C)	1,422 (T=20 °C)		102
	2,3,5,6-Tetrafluoroisonicotinonitrile	4-cyano-2,3,5,6-Tetrafluoropyridine 2,3,5,6-tetrafluoropyridine-4-carbonitrile	2,3,5,6-Тетрафторизоникотонитрил	C <sub>6</sub> F <sub>4</sub> N <sub>2</sub>	16297-07-7	176,073012	99			66-68	
	2,3,5,6-Tetrafluoro-4-methylpyridine	2,3,5,6-Tetrafluoro-4-picoline	2,3,5,6-Тетрафтор-4-метилпиридин	C <sub>6</sub> H <sub>3</sub> F <sub>4</sub> N	16297-14-6	165,090132	99	1,44	1,414		130-131
	3,5,6-Trifluoro-2-picoline	2-methyl-3,5,6-trifluoropyridine	3,5,6-Трифтор-2-пиколин	C <sub>6</sub> H <sub>3</sub> F <sub>3</sub> N			98				
	3,5,6-Trifluoro-2-pyridinamine	2-Amino-3,5,6-trifluoropyridine	3,5,6-Трифтор-2-пиридинамин	C <sub>5</sub> H <sub>3</sub> F <sub>3</sub> N <sub>2</sub>	3534-50-7	148,087429	99			100-103	
	2,3,5,6-Tetrafluoro-4-pyridineamine	4-Amino-2,3,5,6-tetrafluoropyridine	2,3,5,6-тетрафтор-4-пиридинамин	C <sub>5</sub> H <sub>2</sub> F <sub>4</sub> N <sub>2</sub>	1682-20-8	166,077892	98			85-87	
	2,3,6-Trifluoro-4-pyridineamine	4-Amino-2,3,6-trifluoropyridine	2,3,6-Трифтор-4-пиридинамин	C <sub>5</sub> H <sub>3</sub> F <sub>3</sub> N <sub>2</sub>	63489-55-4	148,087429	98	1,5 (lit)	1,497		251,2±35 (lit)
	3,5,6-Trifluoropyridine-2-acetic acid		3,5,6-Трифторпиридин-2-уксусная кислота	C <sub>7</sub> H <sub>4</sub> F <sub>3</sub> NO <sub>2</sub>			98				
	Ethyl (2,3,5,6-tetrafluoro-4-pyridinyl) acetate		Этил (2,3,5,6-тетрафтор-4-пиринидил) ацетат	C <sub>9</sub> H <sub>7</sub> F <sub>4</sub> NO <sub>2</sub>			97				

	Methyl 3,3,3-trifluoro-2-oxopropanoate	3,3,3-trifluoropyruvic acid Methyl ester	Метил 3,3,3-трифтор-2-оксопропионат	C <sub>4</sub> H <sub>3</sub> F <sub>5</sub> O <sub>3</sub>	13089-11-7	156,061229		1,39	1,332-1,334		84-86
	Perfluoropropanoic acid		Перфторпропионовая кислота	C <sub>3</sub> HF <sub>5</sub> O <sub>2</sub>	422-64-0	164,031755		1,561	1,283-1,286		96
CBr <sub>3</sub> F	Tribromo(fluoro)methane		Трибром(фтор)метан	CBr <sub>3</sub> F	353-54-8	270,721403	98	2,765	1,524 – 1,526	-73	106-107
	1,2,3,4,5-Pentafluoro-6-vinylbenzene	Pentafluorostyrene	1,2,3,4,5-Пентафтор-6-винилбензол	C <sub>8</sub> H <sub>3</sub> F <sub>5</sub>	653-34-9	194,103835	98	1,406	1,4445-1,4465		62-63 (p=50 torr)
	Pentafluorobenzenesulphonyl chloride		Пентафторбензолсульфонил хлорид	C <sub>6</sub> F <sub>5</sub> SO <sub>2</sub> Cl	832-53-1	266,569815	98	1,797	1,478-1,48		210-211
	Ethyl (4-Trifluoromethyl-2,3,5,6-tetrafluorophenyl)acetate		Этил (4-Трифторметил-2,3,5,6-тетрафторфенил)ацетат	C <sub>11</sub> H <sub>7</sub> F <sub>7</sub> O <sub>2</sub>	32304-29-3		98				
	2-Fluoro-4-(trifluoromethyl)benzotrile	4-Cyano-3-fluorobenzotrifluoride	2-Фтор-4-(трифторметил)бензонитрил	C <sub>8</sub> H <sub>3</sub> F <sub>4</sub> N	146070-34-0	189,112132	99	1,358 (T=20°C)	1,446 (T=20°C)		178-180
	Pentafluorobenzenethiol	Pentafluorothiophenol	Пентафторбензотиол	C <sub>6</sub> HF <sub>5</sub> S	771-62-0	200,125955	98	1,5	1,463-1,465	-24	143
	1,1'-Sulfanediylbis(pentafluorobenzene)	Bis(pentafluorophenyl)sulphide	1,1'-Сульфанилбис(пентафторбензол)	C <sub>12</sub> F <sub>10</sub> S	1043-50-1	366,17603	98			86-88	124-125
	1-(2-Brom-1,1,2,2-tetrafluoroethoxy)-1,2,2-trifluoroethene	2-Bromotetrafluoroethyl trifluorovinyl ether	1-(2-Бром-1,1,2,2-тетрафторэтокси)-1,2,2-трифторэтан	C <sub>4</sub> BrF <sub>7</sub> O	85737-06-0	276,936221	99				55
	Methylsilane	Monomethylsilane	Метилсилан	CH <sub>3</sub> SiH <sub>3</sub>	992-94-9	46,14414	99,9			156,5	